

ENRICHMENT MEDIA

THIOGLYCOLLATE

General use enrichment medium used for the isolation and culture of large variety of aerobes, anaerobes and microaerophiles. It is used as an enrichment medium for clinical samples. For single use. Supplied in a 5ml tube containing 2 ml of medium.

SELENITE

Selective enrichment broth for *Salmonella spp.* and some *Shigella* species from feces, urine, food and/or water. Inhibits intestinal microbiota proliferation. It is used as an enrichment medium for clinical samples and for industrial use. For single use. Supplied in a 5ml tube containing 2 ml of medium. Occasional presence of light color variations and/or an orange precipitate does not affect the medium viability.

STORAGE CONDITIONS

It is advised to store the enrichment media in a cool, dry place, protected from light, at temperatures between 2-30 °C. Store the tubes containing samples in a cool dark place.

PRECAUTIONS

Sterile product dosed in aseptic conditions. Not suitable for any other application than its intended use. Do not use if there is any sign of packaging breakage or damage. Do not use if any change is observed in the media. Due to the intended use of the product, the end users belong to the Health sector, therefore the used product should be discarded and managed as biological waste using the specific biological waste containers. The product, once used, must be decontaminated or treated as contaminated medical waste. It is the responsibility of each laboratory to handle, treat and dispose of waste according to regulations.

Symbol glossary:

QTY Quantity	CE CE marking	Keep away from sunlight	LOT Batch code	Storage temperature
REF Catalogue number	Do not use if Package is damaged	Manufacturer	Do not re-use	STERILE A Sterile by Aseptic Technique
Use-by-date	Consult instructions for use on the website www.deltalab.es/eifus	IVD In Vitro Diagnostic Medical Device		

In case of a serious incident* related to the product, notify to Deltalab, S.L. as well as the competent authority of the State in which the user is established. *A "serious incident" is understood as one that entails the death, or serious deterioration of the health of the patient or user or a serious threat to public health.