

BED PAN AND URINE COLLECTION BOTTLE

**For the collection of human biological fluids (urine and feces). Clean and disinfect before use. No sterilised.
The product is not manufactured with natural rubber latex or dried natural rubber.
Medical device.**












Intended use:

The bed pan is intended for the collection of urine and/or feces in bedridden patients or with reduced mobility. The urine collection bottle is intended for urine collection in bedridden male patients. Both the pan and the bottle are Intended for hospital and/or home use.

General precautions:

1. Not suitable for any other application than its intended use.
2. The use of the pan and bottle presents no difficulty for the intended end users.
3. Do not use the pan or bottle if the package is damaged.
4. It should be cleaned with detergent and disinfected with bleach or other disinfectants, before each use. The pan is suitable for use in disinfectant washers.
5. The life of the pan is determined by wear and deterioration due to use. Discard when colour changes, roughness, deformations, etc. are detected.
6. Store away from heat sources.

Symbol glossary:

 REF	Catalogue number	 LOT	Batch code		CE marking
 MD	Medical Device	 QTY	Quantity		Use-by-date
	Keep away from sunlight		Manufacturer		Do not use if Package is damaged
 UDI	Unique device identifier		Consult instructions for use on the website www.deltalab.es		

In case of a serious incident* related to the product, notify to Deltalab, S.L. as well as the competent authority of the State in which the user is established. *A "serious incident" is understood as one that entails the death, or serious deterioration of the health of the patient or user or a serious threat to public health.